






# BACK BELT

## THE FIBER

Xenia garments and devices are made with technologically advanced fabrics and yarns, with the addition of ceramic powders capable of reflecting the far infrared, also called fir. The far infrared frequencies are biocompatible, i.e. Compatible with all living organisms and induce a resonance effect with the water known as the "heat effect", thanks to which microcirculation is optimised.


Recommended for sportsmen, for those who drive many hours a day or for those who carry out jobs that require heavy physical effort, thanks to the 3D processes with a massaging effect, the lumbar belt promotes the reactivation of circulation in the lumbar-sacral region. Furthermore, the transpiration zones support the correct thermoregulation. It offers postural support and gives greater help to the affected area. Lastly, the lumbar belt relieves pain caused by inflammation and facilitates recovery in the re-education and post-operative rehabilitation phase. In all the aforementioned cases, the use of Xenia garments does not exclude or replace what is medically prescribed.



F  
I  
R

<b>Code:</b>	<b>Size</b>	<b>Colour: BLACK</b>
YCL0101	XS/S	
YCL0102	M/L	
YCL0103	XL/XXL	
YCL0104	3XL	

## COMPOSITION

92% PA, 8% EA.  
100% of PA is  byon

## TIPS FOR WASHING AND MAINTENANCE FOR BLACK GARMENTS

At low temperatures, both by hand and in the washing machine.  
It is advisable to: wash in max 40° water; wash separately

[WWW.XENIAGLOBAL.IT](http://WWW.XENIAGLOBAL.IT)



Chirurgia  
della  
Mano  
Savona



100% MADE IN ITALY

**Xenia**  
global

**TESSILE**  
BENESSERE

INFORMATION MATERIAL IS INTENDED EXCLUSIVELY TO SPECIALIZED PERSONNEL.  
REPRODUCTION AND DIFFUSION PROHIBITED

Copyright © 2023 All rights reserved - Xenia srl - Piazza della Repubblica, 28 - 20124 Milan Italy

IMPRESA  
ASSOCIATA  
**avedisco**